



Research Phase Contrast Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A19.2602
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Research
- Certification: CE, Rohs
- Drawtube: Binocular
- Theory: Phase Contrast Microscope
- Nosepiece: Quadruple
- Light Source: Transmitting Light, Kohler Illumination System: 12V/20W Halogen
- Highlight: **phase contrast microscopes, inverted phase contrast microscope**



for more products please visit us on cnoec.com

Product Description


Research Inverted Phase Contrast Microscope With Coaxial Coarse A19.2602

Specifications:

A19.2602 Inverted Phase Contrast Microscope Specification		
Head	Seidentopf Trinocular Head, 45° Inclined, Interpupillary Distance 50-75mm, Light Distribution (binocular tube: photo tube): Both 100:0 and 20:80	
Eyepiece	WF10x/20mm	
	WF16x/16mm	
	10x Reticle Eyepiece (0.1mm)	
	Telescope	
Long Working Distance Objective	Plan Objective LPL4x/0.10 WD=17.5mm	
	Plan Objective LPL25x/0.40 WD=5.0mm	
	Plan Objective LPL40x/0.65 WD=3.0mm	
	Plan Phase Contrast Objective LPL Ph+ 10x/0.25 WD=7.9mm	
	Plan Phase Contrast Objective LPL Ph+ 25x/0.40 WD=5.0mm	
Focus System	Low positioned, coaxial coarse & fine focusing system, coarse 42mm/rot, fine 0.2mm/rot, minimum reading 2µm	
	Quadruple	
Stage	Mechanical Stage 206(X) x 170(Y)mm, with attached mechanical stage moving range 128(X) x 80(Y)mm	
Culture Dish Holder	65mm	
	35mm	
Condenser	Turnable Relief Phase Contrast Condenser N.A.0.35 W.D. =45mm	
Phase Contrast	Phase Contrast Slider 10x	
	Phase Contrast Slider 25x	
	Phase Contrast Slider 40x	
Filter	Blue, Green & Gray	
Light Source	Transmitting light, Kohler Illumination System: 12V/20W Halogen	
Adapter	C-Mount Adapter 1x	
Power Supply	85-265V, 50/60Hz	
Optional Accessories		
	Item No.	
Long Working Distance Objective	Plan Objective LPL2.5x/0.07, WD=29.7mm	A52.2608-2.5
	Plan Objective LPL4x/0.10, WD=17.5mm	A52.2608-4
	Plan Phase Contrast Objective LPL Ph+ 10x/0.25 W.D.=7.9 mm	A5C.2630-10
	Plan Phase Contrast Objective LPL Ph+ 40x/0.25 W.D.=3.0 mm	A5C.2630-40
Culture Dish Holder	54mm	A54.2606-54
Stage	Attached Mechanical Stage for 96 holes Plate	A54.2601
Adapter	C-Mount Adapter 0.6x (Focus Adjustable)	A55.2601-B06
Digital SLR Camera Adapter	For Nikon	A55.2622-N
	For Sony	A55.2622-S
	For Canon	A55.2622-C
Light Source	Upgrade to 6V30W Halogen	A56.2650-30W
Other	Micrometer 0.1mm	A58.2601-M



www.cnoec.com



www.cnoec.com



Focusing Hand Wheel

Quadruple Nosepiece



Attached Mechanical Stage for 95mm Plate

C-Mount Adapter 2-Inch Adjustable

C-Mount Adapter 1.25-Inch



65mm Culture Dish Holder

60mm Culture Dish Holder

35mm Culture Dish Holder



Turnable Relief Phase Contrast Condenser

Slide Relief Phase Contrast Condenser

Kohler Illumination System



Phase Contrast Objectives

Ordinary Phase Contrast Objective

Relief Phase Contrast Objective



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China